

Style SEMITRANSSPARENT



semitransparent

Fine: elegant, semitransparent appearance



elegant

Stylish: feminine esthetics; also available in four trendy colors



perfect shaping

Attractive: shapes leg and buttocks optimally




Matchcode SSE		A-D Calf		A-G Thigh	A-T Pantyhose	A-T Maternity
		Soft top	Grip top	Grip top	Soft top	Soft top
CCL	1	●		●	●	●
	2	●	●	●	●	●
Toe		●	●	●	●	●
		●	●	●	●	●
Grip top	Sensinova grip top			●		
	Quadrinova knobbed grip top		●	●		
Colors	● black ● caramel ● savanna	●	●	●	●	●
	● skin	●		●	●	●
No stock item, produced on order:	● grey ● marine					
	● true blue* ● cool grey*	●		●	●	●
	● deep lavender* ● happy yellow*					

* Available as of 01.03.2021

MAGIC

		Circumference cm											
		XXS	XXS Plus	XS	XS Plus	S	S Plus	M	M Plus	L	L Plus	XL	XL Plus
CT		60-75	65-80	65-80	75-90	75-90	85-100	85-100	95-110	95-110	105-120	105-120	115-130
CH		75-90	80-95	80-95	85-100	85-100	95-110	95-110	105-120	105-120	115-130	115-130	125-140
CG	without grip top	42-54	47-59	47-59	52-64	52-64	56-68	56-68	61-74	61-74	66-79	66-79	73-85
	with grip top	42-52	47-56	47-56	52-61	52-61	56-65	56-65	61-71	61-71	66-76	66-76	73-85
CF		38-49	41-52	41-52	44-55	44-55	49-60	49-60	53-64	53-64	59-70	59-70	63-78
CE		29-36	33-40	33-40	37-44	37-44	39-46	39-46	41-48	41-48	45-52	45-52	46-53
CD		23-30	27-34	27-34	31-38	31-38	34-41	34-41	37-44	37-44	41-48	41-48	42-50
CC		25-32	29-36	29-36	33-40	33-40	35-42	35-42	38-45	38-45	42-49	42-49	43-51
CB₁		22-26	22-26	24-29	24-29	26-32	26-32	29-35	29-35	31-37	31-37	33-39	33-39
CB		18-20	18-20	20-22	20-22	22-24	22-24	24-26	24-26	26-29	26-29	29-32	29-32
CY		27-31	27-31	28-33	28-33	29-35	29-35	31-37	31-37	32-38	32-38	33-39	33-39
CA		18-20	18-20	20-22	20-22	22-24	22-24	24-26	24-26	26-29	26-29	29-32	29-32

CT



Length cm					
		short	56 - 64		
normal	35 - 40	normal	65 - 74	normal	65 - 74
long	41 - 45	long	75 - 84	long	75 - 84



A-D = Calf; **A-G** = Thigh; **A-T** = Pantyhose; **A-T** Maternity = for pregnant and oversized women:
circumferential measurements cH and cT are replaced by a special adapting trouser component

Export • SIGVARIS AG • Gröblistrasse 8 • 9014 St. Gallen • Phone: +41 71 272 40 00 • Fax: +41 71 272 40 01

Great Britain • SIGVARIS BRITAIN LTD. • 1 Imperial Court • Magellan Close • Andover, SP10 5NT • Phone: +44 1264 326 666 • Fax: +44 1264 369 950

Australia / New Zealand • SIGVARIS (Australia) Pty. Ltd. • 207 Victoria Street • West Melbourne VIC 3003 • Phone: 1300 SIGVARIS (1300 744 827) • Fax: +61 39329 3773 • info.au@sigvaris.com

Middle East • SIGVARIS AG (DMCC Branch) • Unit No: AG-05-C • Silver (AG) Tower • Jumeirah Lakes Towers • Dubai - United Arab Emirates • Phone: +971 4 558 6265 • Fax: +971 4 558 6269 • middle-east@sigvaris.com

www.sigvaris.group